

Receptor Membrane Products for Screening

Caliper offers GPCR, ion channel, nuclear receptor, and other receptor membrane reagents for screening. All receptor membranes are derived from either human recombinant or human cell lines expressing the receptor.

Our ready-to-use preparations come with their own pharmacological evaluation including K_d/B_{max} values and graphical representations of reference curve(s) with K_i values for each lot number. Assay specification sheets outlining basic assay parameters and conditions are also included with each shipment. In addition, these same membrane preparations are used at Caliper's Discovery Alliances & Services division, *in vitro* assay services, to ensure continued QA/QC validation and the highest quality materials. Each Receptor Membrane is available in either 100 or 400 unit sizes.

Receptor Membrane Products

- Adenosine Transporter Receptor Membranes
- Adrenergic, Alpha 2A Receptor Membranes
- Adrenergic, Alpha 2B Receptor Membranes
- Adrenergic, Beta 1 Receptor Membranes
- Adrenergic, Beta 2 Receptor Membranes
- Adrenergic, Beta 3 Receptor Membranes
- Anaphylatoxin, Complement C5a Receptor Membranes
- Androgen (Testosterone) Receptor Membranes
- Angiotensin II, Type I Receptor Membranes
- Benzodiazepine, peripheral Receptor Membranes
- Dopamine, D1 Receptor Membranes
- Endothelin, ETA Receptor Membranes
- Estrogen Receptor Membranes
- Neuropeptide, NPY2 Receptor Membranes
- Nicotinic, aBnTx insensitive Receptor Membranes
- Progesterone
- Serotonin Transporter Receptor Membranes
- Serotonin, 5HT3 Receptor Membranes



To Place An Order

Go online to:
<https://store.caliperls.com>
or call us at:
1-410-712-4410



Corporate Headquarters
68 Elm Street
Hopkinton, MA 01748-1668
www.caliperLS.com

In Vivo Discovery Alliances
5 Cedar Brook Drive
Cranbury, New Jersey 08512 USA
Tel: 1.609.860.0806

In Vitro Services
7170 Standard Drive
Hanover, Maryland 21076-1334 USA
Tel: 1.410.712.4410 / 1.800.543.4141